

Users Guide

Kit Contents:

Reagent List

- A: Exposure Medium Concentrate 2.5 ml (1 unit)
 - B: D-Glucose 1.5 ml (1 unit)
 - C: Bromocresol Purple 1.5 ml (1 unit)
 - D: D-Biotin 2 ml(1 unit)
 - E: L-Histidine 50 µl (1 unit)
 - F: Sterile Distilled Water 80 ml (1 unit)
 - G: Growth Media 10 ml (1 unit)
 - H: 10X Reversion Solution 4.5 ml (1 unit)
 - I: **Reagent I** 50 µl (1 unit)
 - N: 1X Exposure solution 3.5 ml (1 unit)
- Basic Ames ISO strain: TA 100 *Salmonella typhimurium* (1 unit)
4NQO (5 ng/µl), 50 µl (1 unit)

Required Equipment

- Micropipette using disposable sterile tips in the range of 5 to 200 µl.
- An eight-channel multi-pipette (200 µl).
- Pipette Aid (Pipette gun) using 10 ml sterile pipettes tips.
- A 37°C Incubator.
- Parafilm piece (3 units)




Disposable list

- 24-well plate (2 units)
- 384 well plate (3 units)
- 5 ml tube (2 unit)
- 15 ml conical tube (1 unit)
- A biohazard disposable bag

Applications:

- Testing of industrial effluents for presence of possible mutagenic compounds.
- Screening of municipal discharges for spill contamination, improper chemical disposal.
- Routine monitoring of waste water effluent for quality and mutagenicity.
- Screening of recycled potable water supplies for presence of priority pollutants and genotoxins.
- Screening air particulate mater (PM) for sub chronic human health effects.
- Evaluating water and soil samples for elevated levels of personal care product (PCP) residues.





Standard Ames Strains





Strain	Type of Reversion mutation	Standard mutagen	Comment
<i>Salmonella typhimurium</i>: Histidine dependent			
TA97a	Frameshift	9-AA or ICR191	Derived from TA1537 and includes plasmid pKM 101 which induces error-prone DNA repair enzymes to increase sensitivity
TA98	Frameshift	2NF	Derived from TA1538 and includes plasmid pKM 101 which induces error-prone DNA repair enzymes to increase sensitivity
TA100	Base-pair substitution, oxidative	4NQO	Derived from TA1535 and includes plasmid pKM 101 which induces error-prone DNA repair enzymes to increase sensitivity
TA1535	Base-pair substitution, oxidative	4NQO	uvrB repair deficient, rfa mutation increases permeability to mutagens. Sensitive to 3 unique mutagens compared to TA100 (acetaldehyde oxime, 6-mercaptopurine and 1,3-butadiene)
<i>Escherichia coli</i>: Tryptophan dependent			
E. coli WP2 uvrA (pKM101)	Base- pair substitution	4NQO	uvrA deletion mutation eliminates accurate excision repair mechanism, contains plasmid pKM 101 which induces error-prone DNA repair enzymes and increases sensitivity



AMES - 384 ISO™

384 Well Format






Ames-384 ISO™ Advantages



- Reagents, bacteria and other consumable components are supplied ready-to-use in a non-specialized laboratory
- Assay endpoints are easy to read colorimetric changes that require no specialized training
- Procedure uses less reagents, plastics and incubation space compared to traditional Ames tests. More samples can be run with less cleanup.
- Less sample is required for testing. Small sample volumes still produce great results.
- Pre-exposure and bacterial dilution steps ensure mutagen exposure occurs during bacterial log growth phase to improve results.

LEGEND

- 48 hr ● Non-revertant well
- 37 °C ● Revertant well



Sample 384 ISO result plate




Bioactivation (S9 kit add-on)

- Many mutagens must first be metabolized into their reactive form by enzymes
- Depending on the compound under study, bioactivation may be required for detection
- EBPI offers traditional methods of pro mutagen activation through the addition of S9 liver fraction
- EBPI employs a commonly used metabolic activation system which includes post-mitochondrial liver fractions isolated from Sprague Dawley rats, supplemented with cofactors
- The rats are pre-treated with Aroclor 1254 to stimulate enzyme production prior to liver extraction
- This option is offered with all mutagenicity and genotoxicity testing kits for a supplemental cost

S9 Kit Contents

- S9A:** MgCl₂ + KCl solution 0.4 mL (1 unit)
- S9B:** Glucose-6-phosphate solution 80 μL (1 unit)
- S9C:** NADP solution 250 μL (1 unit)
- S9D:** Phosphate buffer 1.5 mL (1 unit)
- S9F:** Lyophilized S9 fraction






Ames-384 ISO™ Procedure




Experiment overview

- Negative control
- Positive control
- Sterility control
- No dilution
- First dilution
- Second dilution
- Third dilution
- Fourth dilution
- Fifth dilution





Exposure Plate
(24 well plate)

0.5 ml



1 X Reversion Plate
(24 well plate)

0.05 ml



3 X (384 well plate)

Note: Prior to using our test kits, we highly recommend the development of individual outlines that are representative of the respective experiment. This outline is only provided as a guideline for one possible method


The image above is an example of a test that utilizes one sample with five distinct dilutions in a triplicate manner.



Detailed Procedure

1. Overnight inoculation and initial bacterial growth

- Always use aseptic techniques for all steps in this procedure
- Add **Reagent V** to nutrient broth, mix and immediately transfer to bacterial bottle
- Shake to dissolve and place in incubator overnight at 37 °C with shaking (if possible)
 - aeration will aid bacterial growth.



2. Bacteria dilution

- The next morning observe bacterial vial for turbidity
- If turbidity is seen proceed with initial OD₆₀₀ measurement
- If turbidity is not seen continue growth in incubator
- Perform bacterial dilution upon OD₆₀₀ measurement to working dilution OD₆₀₀ value given for each strain (eg. TA100 = 0.07)
 - Bacterial dilution encourages log growth phase and increases uptake of solution components

Dilution Calculation

X = Overnight Calculation,

Y = Working Concentration,

Final volume = 3 mL

Volume of overnight bacteria (mL) for dilution = $(Y/X) \times 3 \text{ mL}$


Example:

TA100 was grown overnight with an overnight **OD₆₀₀ = 1.6**

Working dilution of TA100 is **OD₆₀₀ = 0.07**

Volume of overnight TA100 for dilution = $(0.07/1.6) \times 3 \text{ mL} = \mathbf{0.131 \text{ mL}}$

Therefore, you will add **0.131 mL of overnight TA100 bacteria** in **2.869 mL of 1X exposure buffer** to give a final volume of 3 mL.




3. Exposure media preparation

- Prepare exposure media according to guidelines below
- Ensure final solution is well mixed before continuing to next step

Master Exposure Mix (2.5 ml per experiment)


● A) Exposure medium concentrate	2.075 ml
● D) D- Biotin	0.150 ml
● E) L-Histidine	0.025 ml
● B) 40% D-Glucose	+ 0.250 ml
Total	2.50 ml



4. Transfer exposure media to 24 well plate




100 μ l



Master Exposure Mix

Exposure Plate


5. Transfer positive control mutagen to respective wells



- Pipette **50 µl of 4NQO** into assigned positive control wells
- Two positive control wells are used to increase statistical significance of response
- Mix wells by repeatedly reinjecting solution through pipette

6. Transfer undiluted sample to respective wells

- Direct add 800 μ L of undiluted sample to three wells. This will serve as the most concentrated sample of the series
- One well is left without sample or bacteria to check the sterility of the kit reagents (dark blue)




- (Blue) Reagent blank (no sample no bacteria)
- (Black) Undiluted sample wells

Pipette 800 μ L of undiluted samples into assigned wells

7. Transfer diluted samples to respective wells


- There are enough wells in the exposure plate for 6 different concentrations. (1 undiluted + 5 dilutions)
- Each concentration is run in triplicate to enhance the statistical relevance of the results
- Dilution ratios can be changed depending on purpose of assay
 - Initial screening experiments may warrant wider ranges to find effective concentrations



Pipette diluted samples ($X \mu\text{l} \leq 800 \mu\text{l}$) into assigned wells

8. Transfer sterile water to 24-well plate


- Sterile water is added to bring the **FINAL** volume of all wells to 900 μL
- Depending on sample concentrations amount of water added to the sample wells may require adjustment



- 800 μL of sterile water
- 750 μL of sterile water
- 900 μL of sterile water
- Volume of sterile water = 800 μL – **volume of samples**
- 0 μL of sterile water

9. Transfer bacteria to 24-well plate

- Bacteria is added to all wells except the reagent control (sterility control well, denoted by purple)
- **FINAL** volume of all wells after bacterial addition should be 1 mL





Pipette **100 µl of diluted bacteria** into assigned wells

10. Exposure Incubation

- Mix exposure plate prior to incubation or during incubation (if possible)



- Negative control
- Positive control
- Sterility control
- No dilution
- First dilution
- Second dilution
- Third dilution
- Fourth dilution
- Fifth dilution



Incubate at 37°C
For 100 minute

11. Preparation of Reversion solution and Plate

- Near the end of the incubation period, prepare master reversion mix solution as outlined below
- Dispense 300 µL of reversion mix and 2.2 mL of sterile water to all wells of a second 24-well plate




● B) 40% D-Glucose	0.77 ml
● C) Bromocresol purple	1.2 ml
● D) D-Biotin	1.58 ml
● H) 10X Reversion Solution I	3.95 ml
Total	7.5 ml

1X Reversion Plate

Add 0.3 mL of Master reversion mix and 2.2 mL of sterile water into each well of the reversion plate

12. Transfer bacteria from exposure plate to reversion plate


- Negative control
- No dilution
- Third dilution
- Positive control
- First dilution
- Fourth dilution
- Sterility control
- Second dilution
- Fifth dilution



After 100 min incubation, transfer 0.5 mL of exposed bacterial solution from exposure plate into respective wells in the Reversion Plate.


13. Transfer bacteria-reversion mix into 384-well plate

- Transferring 50 µL of bacteria-reversion mix using a multichannel pipette
- 384-well plate is divided into rows of 24 wells
- Each sample is pipetted into two rows (48 wells) to observe reversion
- Each plate can accommodate 8 samples



14. Incubation of 384 well plate

- Incubate 384-well plates at 37 °C for 2 days
- After two days, remove plates and observe the number of wells that have changed from purple to yellow
- Ensure that negative controls have small numbers of reversions, positive controls have many revertants and reagent sterility control does not have any revertants




Incubate at 37°C
For 2 days

Result Interpretation

LEGEND

- 48 hr ● Non-revertant well
- 37 °C ● Revertant well



Sample 384-well result plate

1. Plates were scored visually. Yellow and partial yellow wells are scored as positive. Purple wells are scored as negative.
2. The test is valid if the following criteria are met:
 - a) Observe the 'Blank' (sterility assessment) wells. Proceed only if the blank wells are sterile (purple). If the well is turbid or yellow, the assay may be contaminated or the sample is interfering with the reagents; results will be invalid.
 - b) Average score for negative or background control is ≥ 0 and ≤ 15 revertant wells per 48-well section on day 2.
 - c) Average score for positive (standard mutagen) controls is ≥ 25 revertant wells per 48-well section on day 2.

If any of these criteria are not met, the test is considered invalid!

Result Significance

No. Wells Positive in Background Plate	No. Wells Positive in Treatment Plate			No. Wells Positive in Background Plate	No. Wells Positive in Treatment Plate		
	0.05	0.01	0.001		0.05	0.01	0.001
0	3	6	9	24	33	36	40
1	5	8	12	25	34	37	40
2	7	10	14	26	35	38	41
3	9	12	15	27	36	39	42
4	10	13	17	28	36	40	43
5	12	15	19	29	37	40	43
6	13	16	20	30	38	41	44
7	14	18	21	31	39	42	45
8	16	19	23	32	40	43	46
9	17	20	24	33	41	44	46
10	18	21	25	34	42	44	47
11	19	23	27	35	43	45	47
12	20	24	28	36	43	46	48
13	21	25	29	37	44	46	49
14	22	26	30	38	45	47	49
15	24	27	31	39	46	48	50
16	25	28	32	40	46	48	50
17	26	29	33	41	47	49	50

- Depending of amounts of revertant wells in tested samples and revertant wells in negative controls, different levels of significance can be assigned
- Use the quick reference chart included with your procedure (left) or more advanced statistical methods to assign significance to mutagenicity results

In general, for mutagenicity to be suggested, there must be greater than double the amount of revertant wells compared to controls and evidence of a concentration-dependant dose-response relationship. (lower doses produce lower amounts of revertants)


AMES 384 ISO EXCEL Spreadsheet

The screenshot displays the AMES 384 ISO Excel spreadsheet interface, which is used for data entry, analysis, and visualization of results. The interface is divided into several sections:

- Top Left:** A Paint application window showing a logo for ebpi environmental bio-detection products inc. with a legend for significance levels: 95% significance (blue), 99.9% significance (green), Possibility cytotoxic (less than average of backgrounds) (yellow), and Value is out of range (red).
- Top Right:** A Microsoft Excel window showing the 'Muta-ChromoPlate 96 template' with a 'Plate Entry' dialog box. The dialog box allows for selecting plate descriptions and bacterial strains for each well.
- Center:** A Microsoft Excel window showing the main data entry spreadsheet. The spreadsheet has columns for Experiment ID, Plate Number, Plate Description, Bacterial Strain, Day, and Number of Positively Scored Wells. A 'UseForm1' dialog box is overlaid on the spreadsheet, asking 'How many plates for this experiment?' with a value of 24 entered.
- Bottom:** Two line graphs showing 'Positive Wells' over 'Day' for two experiments: 'Positive Wells Experiment 1 TA 98 (S9)' and 'Positive Wells Experiment 1 TA 97'. The graphs plot the number of positive wells (0 to 100) against the day (0 to 7). The legend for TA 98 includes Background S9, Test Samples 1-9, and Positive Control 2-AA S9. The legend for TA 97 includes Background, Test Samples 1-9, and Positive Control 2-4A S9.

AMES - 384 ISO™

384 Well Format



Kit Options

- Basic Kits with and without or without S9 Activation
- TA100/TA98 Kits with or without S9 Activation
- OECD 471 Bacterial Based Kits
- 5, 10 Sample Kits
- Reagents Only Kits
- All kits/reagents can be modified to meet your requirements
- Ames Express Strains

